

VOLO Assemble Rack

VOLO 4-way Assemble Racks

- 매우 낮은 온도에서도 사용이 가능합니다. (-90℃)
- Available in standard or interlocking versions
- 일반적인 튜브 사이즈 대부분을 수용(from 0.5ml to 50ml)
- 단단한 조립 구조 형성

Product List

VOLO 4-Way Assemble Racks				
Item	Description	Tube Options	Colors	Packaging
FP-101	4-Way Interlock Assemble 171 x 91 x 54 mm	Side 1: 32 x 1.5-2.0 mL Side 2: 32 x 0.5 mL Side 3: 12 x 5 or 15 mL Side 4: 4 x 50 mL	Fluorescent: Blue, Green, Yellow, Orange, Pink	5units / pack 4packs / case
FP-102	4-Way Assemble 5ml 171 x 91 x 54 mm	Side 1: 18 x 5 mL and 32 x 0.5 mL Side 2: 32 x 1.5-2.0 mL Side 3: 12 x 5 or 15 mL Side 4: 4 x 50 mL	Standard: Blue, Green, Red, Yellow	4units / pack 4packs / case
FP-103	4-Way Assemble 50ml 8 x 50ml Tubes 198 x 114 x 95 mm	Side 1: 8 x 50 mL Side 2: 20 x 20 mm Side 3: 24 x 5 or 15 mL Side 4: 40 x 1.5-2.0 mL	Fluorescent: Blue, Green, Yellow, Orange, Pink	5units / pack 2packs / case
FP-104	4-Way Assemble 25mm 12 x 25mm Tubes 198 x 114 x 95 mm	Side 1: 12 x 25 mm Side 2: 20 x 20 mm Side 3: 40 x 1.5-2.0 mL Side 4: 24 x 5 or 15 mL	Standard: Blue, Green, Red, Orange, Violet	5units / pack 2packs / case



FP-101



FP-102



FP-103

VOLO Assemble Rack and Sliplock

Product List

VOLO 4-Way Assemble Racks with Tray				
Item	Description	Tube Options	Colors	Packaging
FP-201	4-Way Assemble with Tray Standard Colors 243 × 102 × 41 mm	12 x 0.5 mL 9 x 1.5-2.0 mL 6 x 15 mL 2 x 50 mL	Standard: Green, Blue, Red, Yellow	1tray / pack 5packs / case
FP-202	4-Way Assemble with Tray Fluorescent Colors 243 × 102 × 41 mm	12 x 0.5 mL 9 x 1.5-2.0 mL 6 x 15 mL 2 x 50 mL	Fluorescent: Green, Yellow, Orange, Pink	1tray / pack 5packs / case



FP-201



FP-202

VOLO Sliplock Assemble				
Item	Description	Tube Options	Colors	Packaging
FP-301	Sliplock Assemble 4 Racks and 2 Clear Lids	16 x 0.5 mL tubes per rack. 24 x 1.5/2.0 mL tubes per rack. 60 x 0.1/0.2 mL PCR tubes per rack. 12 x 8-strip or 11 x 12-strip PCR tubes per 2 racks. 1 x non-skirted or semi-skirted 96-well PCR plate per 2 racks.	Green, Blue, Red, Yellow	1workstation / pack 4packs / case



FP-301

VOLO Revolver Assemble

Product List

VOLO Revolver Assemble				
Item	Description	Tube Options	Colors	Packaging
FP-302	Revolver Assemble 257 x 108 x 108 mm	Side 1: 6 x 0.5 mL and 12 x 1.5/2.0 mL Side 2: 8 x 5 or 15 mL Side 3: 2 x 50 mL Side 4: 48 x 0.2 mL	Green, Blue, Red, Yellow	4rack blocks / pack 4packs / case



FP-302

VOLO Assemble Rack(32, 64, 96)

Product List

VOLO Assemble Rack 32-Place(4 x 8 Format)				
Item	Description	Tube Options	Colors	Packaging
FP-401	Rack with Clear Lid 169 × 81 × 48 mm	Side A: 32 x 1.5/2.0mL + 8 x 0.5 mL Side B: 32 x 0.5 mL + 8 x 1.5/2.0 mL	Assorted Fluorescent	5racks with lids / pack 4packs / case
FP-402	Rack Only 169 × 79 × 28 mm	Side A: 32 x 1.5/2.0mL + 8 x 0.5 mL Side B: 32 x 0.5 mL + 8 x 1.5/2.0 mL	Assorted Fluorescent	5 racks/unit, 4 units/case
VOLO Assemble Rack 64-Place(8 x 8 Format)				
Item	Description	Tube Options	Colors	Packaging
FP-403	Rack with Clear Lid 132 × 130 × 45 mm	Side A: 64 x 1.5/2.0 mL Side B: 64 x 0.5 mL	Assorted Fluorescent	5racks with lids / pack 4units / case
FP-404	Rack Only 132 × 127 × 31 mm	Side A: 64 x 1.5/2.0 mL Side B: 64 x 0.5 mL	Assorted Fluorescent	5racks / unit, 4units / case
VOLO Assemble Rack 96-Place(8 x 12 Format)				
Item	Description	Tube Options	Colors	Packaging
FP-405	Rack with Clear Lid 213 × 118 × 50 mm	Side A: 96 x 1.5/2.0 mL Side B: 96 x 0.5 mL	Assorted Standard	5racks with lids / pack 2packs / case
FP-406	Clear Lid Only 206 × 118 × 29 mm	Side A: 96 x 1.5/2.0 mL Side B: 96 x 0.5 mL	Assorted Standard	1 lid / pack 10packs / case
FP-407	0 Rack Only 213 × 113 × 27 mm	Side A: 96 x 1.5/2.0 mL Side B: 96 x 0.5 mL	Assorted Standard	5racks / pack 2packs / case



FP-401



FP-402




FP-403



FP-404



FP-405



FP-406



FP-407