


VOLO CryoVial

Innovative bi-injection system

- Onebody Cap: Polyethylene+silicone o-ring
- Use TPE(Thermoplastic elastomer) material
- Provide a 100% leak-proof seal
- No contamination, excellent screwing grip

Vial Picker

- Facilitates the removal of cryotubes from VOLO Cryo boxes and workstations



Locking System

- Locking system with a base of polycarbonate workstation system
- Skirted tubes only



All Volo Cryovials are printed barcode

- Pre-printed barcode on every tube
- Readable both ways, by all scanner brands

White writing surface

- For specimen identification

Certification

- Free of human DNA, RNase, DNase, ATP and PCR inhibitors
- Tubes in medical grade polypropylene
- Screw caps in medical grade polyethylene, silicone- injected
- USP Class VI, non-cytotoxic & non-haemolytic
- Sterile (Beta-radiation, SAL 10⁻⁶)
- Autoclavable (open) at 121 °C for 20 minutes
- CE-IVD marked
- Certified IATA PI650
- Double traceability : 2D DataMatrix is available
- Working temperature: up to - 196 °C in liquid nitrogen in the gas phase
- Manufactured in Clean Room Class ISO 7(UNI EN ISO 14644-1). Class 10000(US FED STD 209E)
- Made in Europe



sterile



USP
Classe VI**

VOLO CryoVial

Product List

Internal Cap							
Item	CR-101	CR-102	CR-103	CR-104	CR-105	CR-106	CR-107
Volume	1ml	2ml	2ml	4ml	4ml	5ml	5ml
Bottom	skirted	round	skirted	round	skirted	round	skirted
Packaging	400 / case(50/pack)						




External Cap						
Item	CR-201	CR-202	CR-203	CR-204	CR-205	CR-206
Volume	1ml	2ml	2ml	3ml	4ml	5ml
Bottom	skirted	round	skirted	skirted	skirted	skirted
Packaging	400 / case(50/pack)					



Product List

Cap Disk						
Item	CR-301	CR-302	CR-303	CR-304	CR-305	CR-306
Color	blue	white	yellow	purple	red	green
Packaging	1000 / pack					



Cryo boxes								
Item	CR-401	CR-402	CR-403	CR-404	CR-411	CR-412	CR-413	CR-414
Size	1,2ml in/Ex	1,2ml in/Ex	3,4,5ml in/Ex	1,2 ml in	1,2ml in/Ex	1,2ml in/Ex	3,4,5ml in/Ex	1,2 ml in
Hole	25	81	81	100	25	81	81	100
Barcode	Non	Non	Non	Non	Barcoded	Barcoded	Barcoded	Barcoded
Packaging	48ea / case	40ea / case	32 / case	40 / case	48ea / case	40ea / case	32 / case	40 / case



Product List

Vial Picker

Item	CR-420
Packaging	600ea / case



Workstation

Item	CR-601	CR-602	CR-603
Color	Blue	Green	Red
Hole	40	40	40
Size	12.5 x 10 x 2.5	12.5 x 10 x 2.5	12.5 x 10 x 2.5
Packaging	1ea / pk	1ea / pk	1ea / pk



2D - Barcode System

