


VOLO Deep Well Plates

- Pure Polypropylene으로 제작되어 샘플의 장기간 보관 시 온전한 샘플 보존을 가능하게 합니다.
- Certified free of RNase, DNase, DNA, PCR inhibitors
- 알파벳 및 숫자가 표기되어 있어 샘플 위치 추적이 용이합니다. (알파벳:A-H, 숫자:1-12)
- Centrifugation to 4,000 x G
- HTS-friendly: 대부분의 liquid-handling 장비와 호환됩니다.


Product List

VOLO Deep Well Plates


Item	Volume	Description	Packaging
DP-101	1.1ml	Deep Well Plate, Round Bottom	5units / pack 10packs/ case
DP-101S	1.1ml	Deep Well Plate, Round Bottom, Sterile	5units / pack 10packs/ case
DP-102	1.1ml	Deep Well Plate, Square, V-Bottom	5units / pack 10packs/ case
DP-103	2.0ml	Deep Well Plate, Round Bottom	5units / pack 10packs/ case
DP-103S	2.0ml	Deep Well Plate, Round Bottom, Sterile	5units / pack 10packs/ case
DP-104	2.0ml	Deep Well Plate, Square, V-Bottom	5units / pack 10packs/ case
DP-104S	2.0ml	Deep Well Plate, Square, V-Bottom, Sterile	5units / pack 10packs/ case
DP-105	240uL	Deep Well Plate, V-Bottom	5units / pack 10packs/ case
DP-105S	240uL	Deep Well Plate, V-Bottom, Sterile	5units / pack 10packs/ case



DP-101



DP-102



DP-103



DP-104



DP-105

VOLO Reservoir

Item	Volume	Description	Packaging
DP-204	20ml	12-Channer Reservoir	5units / pack 4packs/ case



DP-204

VOLO Deep Well Sealing Foil and Mats


VOLO Sealing Foil and Mats

- 실리콘 매트는 재사용이 가능하며, 샘플 추출을 위해 구멍을 뚫을 수 있습니다.
- 넓은 온도 범위 사용 가능: -80℃~120℃
- 장기간 샘플 보관에 용이(DMSO를 포함한 다른 용액에 대해 높은 내구성)



Product List

VOLO Deep Well Sealing Foil and Mats			
Item	Well Fit	Description	Packaging
DP-201	96well Round or Square	VOLO XR Foil for Extreme Conditons (129.0 x 78.3mm)	5units / pack 10packs/ case
DP-202	96 Round Well	Sealing Mat for 96-W Deep Well Plate, Round Wells	5units / pack 10packs/ case
DP-202S	96 Round Well	Sealing Mat for 96-W Deep Well Plate, Round Wells, Sterile	5units / pack 10packs/ case
DP-203	96 Square Well	Sealing Mat for 96-W Deep Well Plate, Square Wells	5units / pack 10packs/ case
DP-203S	96 Square Well	Sealing Mat for 96-W Deep Well Plate, Square Wells, Sterile	5units / pack 10packs/ case



VOLO Deep Well Tubes & Rack Systems


VOLO Deep Well Tubes & Rack Systems

- HTS-friendly: 대부분의 liquid-handling 장비와 호환됩니다.
- SBS footprint: 8 x 12 format
- 알파벳 및 숫자가 표기되어 있어 샘플 위치 추적이 용이합니다. (알파벳:A-H, 숫자:1-12)
- 샘플 보관, 희석, 혼합, 스크리닝 등 다양한 용도로 사용 가능합니다.

Product List

VOLO Deep Well, Non-Racked Tubes

Item	Volume	Description	Packaging
DP-001	1.2ml	Individual Tubes	1000tubes / pack 10 packs / case
DP-002	1.2ml	Strips of 8 Tubes	125 Strips / pack 10packs / case
DP-003	1.2ml	Strips of 12 Tubes	80 strips / pack 10packs / case



DP-001



DP-002



DP-003

VOLO Deep Well, Cpa Strips

Item	Description	Packaging
DP-020	Strips of 8 Caps	125 Strips / pack 10packs / case
DP-020S	Strips of 8 Caps / Sterile	125 Strips / pack 10packs / case
DP-030	Strips of 12 Caps	80 Strips / pack 10packs / case
DP-030S	Strips of 12 Caps / Sterile	80 Strips / pack 10packs / case



DP-020



DP-030

VOLO Deep Well Tubes & Rack Systems

Secure Vessels for Storage & Reaction

- Rack System 제품에는 투명한 뚜껑이 함께 제공됩니다.
- Rack System 제품은 Robotic workstation이나 장기 샘플 보관에 용이합니다.
- Individual tube, Strips of 8 tubes, Strips of 12 tubes 3가지로 제품이 생산됩니다.

Product List

VOLO Deep Well, Racked Tubes with Lid			
Item	Volume	Description	Packaging
DP-301	1.2ml	Individual Tubes, Racked	10Racks / Pack 10Packs / Case
DP-301S	1.2ml	Individual Tubes, Racked / Sterile	10Racks / Pack 10Packs / Case
DP-302	1.2ml	Strips of 8 Tubes, Racked	10Racks / Pack 10Packs / Case
DP-302S	1.2ml	Strips of 8 Tubes, Racked / Sterile	10Racks / Pack 10Packs / Case
DP-303	1.2ml	Strips of 12 Tubes, Racked	10Racks / Pack 10Packs / Case
DP-303S	1.2ml	Strips of 12 Tubes, Racked / Sterile	10Racks / Pack 10Packs / Case



DP-301



DP-302



DP-303