

VOLO Freeze Rack

- 상온에서 사용할 경우 4℃ 미만의 시료 온도를 약 4시간 동안 유지(덮개 사용 시)
- 온도에 따른 색상 변화로 육안으로 샘플의 온도 상태를 확인 할 수 있음
- VOLO Freeze Rack의 색상변화는 7℃에서 변화폭이 가장 큼

Product List

VOLO Freeze PCR Racks

Item	Description	Cold to Warm Color Change	Est. Time Below 4℃	Packaging
VF-005	qPCR Rack	Blue to White	2.5 Hours	2racks / pack 5packs / case
VF-003	PCR Rack	Purple to Pink	3+ Hours	2racks / pack 5packs / case
VF-004	PCR Rack	Green to Yellow	3+ Hours	2racks / pack 5packs / case

Tube and plate capacity

- qPCR: 64 x 0.1 mL Corbett Rotor-Gene tubes, 16 x 0.2 mL tubes or strips,
- PCR: 96 x 0.1 and 0.2 mL tubes and strip tubes; plates up to 96-well

VOLO Freeze PCR, SBS

Item	Description	Cold to Warm Color Change	Est. Time Below 4℃	Packaging
VF-006	PCR, SBS Rack	Purple to Pink	2.5 Hours	2racks / pack 5packs / case
VF-007	PCR, SBS Rack	Green to Yellow	2.5 Hours	2racks / pack 5packs / case

SBS PCR

- SBS footprint for use with liquid handling equipment
- 96-well PCR plates, Strip tubes, 0.2ml tube 모두 사용 가능
- 뚜껑이 없는 제품

VOLO Freeze Racks, MCT


Item	Description	Cold to Warm Color Change	Est. Time Below 4℃	Packaging
VF-001	MCT Rack	Purple to Pink	3.5+ Hours	1rack / pack 5 packs / case
VF-002	MCT Rack	Green to Yellow	3.5+ Hours	1rack / pack 5 packs / case

Tube and plate capacity


- 24 x 0.5, 1.5, 2.0 mL snap or screw-cap style tubes.



VF-005




VF-003




VF-004




VF-006



VF-007



VF-001



VF-002