

VOLO PCR Strip

- Angularly attached cap: Tube와 Cap의 연결 부위를 강화하여 주위 다른 cap, tube와의 간섭 제거
- Ultra thin wall: 0.25mm 두께의 tube wall을 통하여 효율적인 열 전달과 Crush 방지
- Certified free of Dnase, Rnase, DNA, and PCR inhibitors.
- Medical grade의 Polypropylene으로 제작하여 real time PCR 시 불필요한 Background signal을 제거

VOLO Indi-H Strip

- Robust interwell linkages: 강화된 연결 부위를 통해 Strip Tube가 쉽게 휘어지던 문제 해결

VOLO Indi-S Strip

- Customizable strips: 손쉬운 Tube 분리로 맞춤형 실험 가능

VOLO One-shot Strip

- Hinged attached cap: 빠른 PCR 실험과 효율적인 Scale up 가능
- Optically clear cap을 통하여 real-time PCR 최적화

Product List

VOLO Indi-H Strip

Item	Volume	Description	Packaging
PP-101	0.2ml	Strips of 8 / flat cap / optically clear	120 strips / pack 10 packs / case
PP-102	0.1ml	Strips of 8 / flat cap / optically clear	120 strips / pack 10 packs / case



PP-101



PP-102

VOLO Indi-S Strip

Item	Volume	Description	Packaging
PP-103	0.2ml	Strips of 8 / dome cap / clear	120 strips / pack 10 packs / case
PP-104	0.2ml	Strips of 8 / flat cap / optically clear	120 strips / pack 10 packs / case
PP-109	0.1ml	Strips of 8 / flat cap / optically clear	120 strips / pack 10 packs / case



PP-103




PP-104



PP-109

VOLO one-shot Strip

Item	Volume	Description	Packaging
PP-105	0.2ml	Strips of 8 / flat cap / optically clear	120 strips / pack 10 packs / case



PP-105

VOLO PCR Strip

Product List

VOLO Strip Tubes with Caps			
Item	Volume	Description	Packaging
PP-106	0.2ml	Strips of 8 / flat cap / clear	120 strips / pack 10 packs / case
PP-107	0.2ml	Strips of 8 / dome cap / clear	120 strips / pack 10 packs / case
PP-110	0.1ml	Strips of 8 / flat cap / clear	120 strips / pack 10 packs / case
PP-111	0.2ml	Strips of 12 / flat cap / clear	80 strips / pack 10 packs / case
PP-112	0.1ml	Strips of 12 / dome cap / clear	80 strips / pack 10 packs / case



PP-106



PP-107



PP-110



PP-111



PP-112



PP-121



PP-122



PP-123



PP-124



PP-125



PP-126



PP-127

VOLO Strip Tubes			
Item	Volume	Description	Packaging
PP-121	0.2ml	Strips of 8 Tubes / clear	125 strips / pack 10 packs / case
PP-122	0.1ml	Strips of 8 Tubes / clear	125 strips / pack 10 packs / case
PP-123	0.1ml	Strips of 12 Tubes / clear	80 strips / pack 10 packs / case

VOLO Strip Caps			
Item	Caps	Description	Packaging
PP-124	Flat	Strips of 8 Caps / clear	125 strips / pack 10 packs / case
PP-125	Dome	Strips of 8 Caps / clear	125 strips / pack 10 packs / case
PP-126	Flat	Strips of 12 Caps / clear	80 strips / pack 10 packs / case
PP-127	Dome	Strips of 12 Caps / clear	80 strips / pack 10 packs / case

VOLO Individual PCR Tubes

- Ultra thin wall: 0.25mm 두께의 tube wall을 통하여 효율적인 열전달과 Crush 방지
- Certified free of Dnase, Rnase, DNA, and PCR inhibitors.
- 증발에 강한 특성과 열고 닫기 쉬운 Cap을 통해 샘플의 손실을 예방합니다.

Corbett-Type Tubes and Caps

- QIAGEN, Corbett Rotor-Gene 장비와 호환 가능합니다.

Product List

VOLO Individual PCR Tubes			
Item	Volume	Description	Packaging
PP-130	0.2ml	Capless Tube	1000 Tubes / pack 10 packs / case
PP-131	0.2ml	dome cap / clear	1000 Tubes / pack 10 packs / case
PP-132	0.1ml	Frosted flat cap / clear	1000 Tubes / pack 10 packs / case
PP-133	0.2ml	flat cap / clear	1000 Tubes / pack 10 packs / case
PP-134	0.5ml	dome cap / clear	1000 Tubes / pack 10 packs / case
PP-135	0.5ml	Frosted flat cap / clear	1000 Tubes / pack 10 packs / case



Clear cap on
PP-132 is ideal
for labeling.

PP-132

VOLO Corbett-Type Strip Tubes and Caps			
Item	Volume	Description	Packaging
PP-108	0.1ml	4 Tube Strip and 4 Cap Strip / Clear	250 Strips / pack 10 packs / case



PP-108