


VOLOSeal EnduTubes

- Bead-beating의 용도에 맞게 특수 제작된 독점 resin으로 제조
- 정교한 O-ring이 튼튼한 Seling 제공
- Pure polypropylene(FDA 21 CFR 177.1520)
- Graduated 100uL: 샘플 볼륨을 보다 쉽게 확인 가능
- Self-Standing
- Certified free of RNase, DNase, DNA, PCR inhibitors

Product List

VOLO EnduTubes & Caps			
Item	Volume	Description	Packaging
PT-501	2.0ml	Graduated EnduTubes & Caps, Assembled	200units / pack 5 packs / case



VOLOSeal Screw Caps with O-Rings

VOLOSeal Grip Caps

- Pure polypropylene(FDA 21 CFR 177.1520 & USP Class VI)
- 고무 그립으로 제작되어 냉동상태에서도 쉽고 열고 닫을 수 있습니다.

Product List

VOLOSeal Grip Caps		
Item	Description	Packaging
PT-403	Rubber ribs on Grip Caps	500units / pack 10 packs / case



VOLOSeal Standard Screw Caps

VOLOSeal Standard Screw Caps

- Available with and with out attached tether
- 합성 O-ring으로 제조(PT-303, PT-304 제외)

VOLOSeal Standard Screw Caps		
Item	Description	Packaging
PT-301	Standard Caps	500units / pack 10 packs / case
PT-302	Standard Cap with Tether	500units / pack 10 packs / case
PT-303	Insertable Cap with Silicone O-ring	500units / pack 10 packs / case
PT-304	Insertable Tethered Cap with Silicone O-ring	500units / pack 10 packs / case
PT-305	Cap Inserts	500units / pack 10 packs / case

PT-301 PT-302



PT-303 PT-304



PT-305



VOLOSeal Screw MicroTubes

VOLOSeal Screw MicroTubes

- Pure polypropylene(FDA 21 CFR 177.1520 & USP Class VI)
- 넓은 온도 범위 사용 가능: -196℃~121℃ (액체질소탱크 사용 불가)
- Rated to 2,000 x G
- Graduated versions: 100uL 마다 눈금이 있어, 샘플 용량을 쉽게 확인 가능
- Certified free of RNase, DNase, DNA, PCR inhibitors
- Skirted or conical bases
- Ideal for boiling applications

Product List

VOLOSeal Clear Tubes, Graduated

Item	Volume	Description	Packaging
PT-202	1.5ml	Conical Base	500units / pack 10 packs / case
PT-203	2.0ml	Conical Base	500units / pack 10 packs / case
PT-201	2.0ml	Skirted Base	500units / pack 10 packs / case

Graduated



VOLOSeal For Manufacturers: Clear Tubes, Non-Graduated

Item	Volume	Description	Packaging
PT-204	1.5ml	Conical Base	500units / pack 10 packs / case
PT-205	2.0ml	Conical Base	500units / pack 10 packs / case
PT-206	2.0ml	Skirted Base	500units / pack 10 packs / case

Non-graduated



VOLOSeal Screw MicroTubes

Product List

VOLOSeal Clear Tubes, Non-Graduated			
Item	Volume	Description	Packaging
PT-101	0.5ml	Skirted Base, Ribbed	500units / pack 10 packs / case
PT-102	0.5ml	Conical Base	500units / pack 10 packs / case
PT-103	0.5ml	Skirted Base	500units / pack 10 packs / case
PT-104	1.5ml	Conical Base, Ribbed	500units / pack 10 packs / case
PT-105	1.5ml	Skirted Base, Ribbed	500units / pack 10 packs / case
PT-106	1.5ml	Skirted Base	500units / pack 10 packs / case
PT-107	2.0ml	Conical Base, Ribbed	500units / pack 10 packs / case
PT-108	2.0ml	Skirted Base, Ribbed	500units / pack 10 packs / case

Ribbed



VOLOSeal Amber Tubes			
Item	Volume	Description	Packaging
AMT-010	1.5ml	Conical Base, Graduated	500units / pack 10 packs / case
AMT-020	2.0ml	Conical Base, Graduated	500units / pack 10 packs / case
AMT-030	2.0ml	Skirted Base, Graduated	500units / pack 10 packs / case
AMT-001	0.5ml	Skirted Base, Ribbed, Non-Graduated	500units / pack 10 packs / case
AMT-002	0.5ml	Skirted Base, Non-Graduated	500units / pack 10 packs / case
AMT-003	1.5ml	Skirted Base, Ribbed, Non-Graduated	500units / pack 10 packs / case
AMT-004	1.5ml	Skirted Base, Non-Graduated	500units / pack 10 packs / case
AMT-005	2.0ml	Skirted Base, Ribbed, Non-Graduated	500units / pack 10 packs / case

Amber



VOLOSeal Tube & Cap Assemblies, Sterile

- Pure polypropylene(FDA 21 CFR 177.1520 & USP Class VI)
- 넓은 온도 범위 사용 가능: -196℃~121℃ (액체질소탱크 사용 불가)
- Certified free of RNase, DNase, DNA, PCR inhibitors, Pyrogens

Product List

Tubes + Standard Caps (PT-301)
(50 assemblies / pack, 50packs / case)

Graduated	0.5ml Conical Base	0.5ml Skirted Base	1.5ml Conical Base	1.5ml Skirted Base	2.0ml Conical Base	2.0ml Skirted Base
Non-Graduated	PT-102-C	PT-103-C	-	PT-106-C	-	-
Graduated	-	-	PT-202-C	-	PT-203-C	PT-201-C
Ribbed	-	PT-101-C	-	PT-105-C	-	PT-108-C

Tubes + Standard Tethered Caps(PT-302)
(50 assemblies / pack, 50packs / case)

Graduated	0.5ml Conical Base	0.5ml Skirted Base	1.5ml Conical Base	1.5ml Skirted Base	2.0ml Conical Base	2.0ml Skirted Base
Non-Graduated	-	PT-103-TC	-	PT-106-TC	-	-
Graduated	-	-	PT-202-TC	-	PT-203-TC	PT-201-TC
Ribbed	-	-	-	-	-	-

