


Gel Loading Tips



Sorenson is *THE* trusted manufacturer of gel loading tips. Sorenson's MultiFlex™, MultiFlex E and MiniFlex gel loading tips have become the industry standard and have been used globally in genetic sequencing studies.

MultiFlex Tips (1-200 µl):


- Universal fit on research-grade pipettors
- 83 mm length for access to bottom of narrow vessels and gel wells
- 5 µl graduation for volume QuickCheck
- Available in bulk (round) and racked (round and flat)

MiniFlex Tips (0.1-10 µl):

- Universal fit on ultra-micro pipettors
- 45 mm length for access to bottom of narrow vessels and gel wells
- 2 µl graduation for volume QuickCheck
- Available racked (round and flat)

MultiFlex E Tips (0.5-10 µl):

- Fit Eppendorf and ultra-micro pipettors
- 5 µl graduation for volume QuickCheck
- Microcapillary is ideal for viewing small sample volumes
- Available racked (round and flat)



White racks – Round (0.57 mm O.D.) tips

Yellow racks – Flat (0.37 mm O.D.) tips

Blue racks – Flat (0.17 mm O.D.) tips

Gel Loading Tips

	Cat. No.	Vol. Range	Orifice Shape	Orifice O.D.*	Graduations	Packaging	Pack Type**	Quantity
10 µl MiniFlex 	15100	0.1-10 µl	Round	0.57 mm	2 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	15110	0.1-10 µl	Round	0.57 mm	2 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	15120	0.1-10 µl	Flat	0.37 mm	2 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	15130	0.1-10 µl	Flat	0.37 mm	2 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	17350	0.1-10 µl	Flat	0.17 mm	2 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	17360	0.1-10 µl	Flat	0.17 mm	2 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
10 µl MultiFlex E 	13920	0.5-10 µl	Round	0.57 mm	5 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	13930	0.5-10 µl	Round	0.57 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	13890	0.5-10 µl	Flat	0.37 mm	5 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	13900	0.5-10 µl	Flat	0.37 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	17330	0.5-10 µl	Flat	0.17 mm	5 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	17340	0.5-10 µl	Flat	0.17 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
200 µl MultiFlex 	28480	1-200 µl	Round	0.57 mm	5 µl	Bulk, Non-Sterile	Bulk Pack	1000 Tips/Bag, 10 Bags/Case
	13790	1-200 µl	Round	0.57 mm	5 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	13810	1-200 µl	Round	0.57 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	13760	1-200 µl	Flat	0.37 mm	5 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	13770	1-200 µl	Flat	0.37 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	17310	1-200 µl	Flat	0.17 mm	5 µl	Racked, Non-Sterile	Flex Rack	200 Tips/Pack, 4 Packs/Case
	17320	1-200 µl	Flat	0.17 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case

*.37 mm O.D. tips are for use with 0.4 mm gels; .17 mm O.D. tips are for use with 0.2 mm gels

**Gel Loading Tip Packaging: an individual pack is a tip rack of 200 tips